

## KINETIC TOXICITY DETECTION OF ENVIRONMENTAL SAMPLES


### BioTox™ kit (Aboatox)

#### Introduction

Triathler™ luminometer equipped with automatic injector together with freeze-dried BioTox™ reagents provide fast and sensitive method for determination of the inhibitory effect of all water soluble samples, including suspensions of solid samples, on the light emission of *Vibrio fischeri* NRRL-B 11177. This application note describes the kinetic method for detection the toxicity, overruling the need for extra correction of quenching caused by color and turbidity.

#### Principles of the test

The BioTox™ reagent contains naturally luminescent *Vibrio fischeri* bacteria, producing luciferase as a part of its metabolic pathway. Luciferase specifically reacts with FMNH<sub>2</sub> and catalyzes the following reaction.



Substances affecting any part of the metabolic pathway of the bacteria, directly affect to the amount of light emitted.

#### Assay Procedure on Triathler™

The inhibitions (INH%) caused by the sample dilutions to *V. fischeri* bacteria are determined and EC<sub>50</sub> –values (sample concentration causing 50% inhibition) calculated using linear regression analysis. The sample dilutions and control sample (2% NaCl) are pipetted in the test tubes (500 µL each). The Triathler injector is filled with the *V. fischeri* reagent. The first tube containing the control sample (2% NaCl) is placed in Triathler and the measurement and injection of 500 µL of the reagent

initiated simultaneously. The measurement is continued on kinetic mode for 5 seconds. The tube is placed in incubator at 15°C and the same procedure repeated for each sample. After incubation time of 15 minutes the end-point readings from the control and each sample are being measured.

#### Procedural Notes:

The tubes shall be shaken before the measurements in order to obtain homogenous suspension. Any incubation time 5 – 30 minutes can be used instead of 15 minutes.

#### Results

Inhibitions (INH%) of each sample dilution are calculated according to the equations below and plotted on log-log scale. The EC<sub>50</sub>-value is determined by using standard linear regression analysis.

$$\text{KF} = \frac{\text{IC}_{15}}{\text{IC}_0} \qquad \text{INH\%} = 100 - \frac{\text{IT}_{15}}{\text{KF} \times \text{IT}_0} \times 100$$

Where:

KF = Correction factor.

IC<sub>15</sub> = Luminescence intensity of the control after contact time (15 min) in CPS.

IC<sub>0</sub> = Maximum CPS value of the control during the 5 second kinetic measurement.

IT<sub>15</sub> = Luminescence intensity of test sample after contact time (15 min) in CPS.

IT<sub>0</sub> = Maximum CPS value of the sample during the 5 second kinetic measurement.

#### Conclusions

The BioTox kit together with Triathler luminometer and automatic injector offer the fastest and most easy-to-use method for detection of toxicity of colorful and turbid samples, including water suspensions of solids (soil, waste and sediment). Also clear water samples can be measured.

#### Materials required

Triathler™ luminometer

BioTox™ kit (Code 1243-500)

Luminometer cuvettes or similar test tubes.

Incubator to maintain the reagent/samples at 15°C.

### Product Information

**Finland and International** Hidex Oy tel +358-2-2750557 fax +358-2-2410075 [info@hidex.com](mailto:info@hidex.com) **Australia** Skudtek Scientific Pty Ltd Tel +613 9793 2929 [skudtek@ozemail.com.au](mailto:skudtek@ozemail.com.au) **Austria** HVD Life Science Vertriebs GmbH Tel +43 1 982 9526 [office1@hvdgmbh.com](mailto:office1@hvdgmbh.com) **France** ScienceTec Tel +3316446 2400 [info@sciencetec.fr](mailto:info@sciencetec.fr) **Germany** Zinsser Analytic GmbH Tel +49 69 789 1060 [info@zinsser-analytic.com](mailto:info@zinsser-analytic.com) **Japan** Sowa Trading Co Inc. Tel +81 3 3862 2700 [tokyo@sowa-trading.co.jp](mailto:tokyo@sowa-trading.co.jp) **Korea** SIMS Co.Ltd Tel +82 2 576 9727 [dmllee@sims.com](mailto:dmllee@sims.com) **Middle East** HVD Vertriebs GmbH tel +30 10 960 06 87 [office@hvd.gr](mailto:office@hvd.gr) **Russia** Pribori Oy c/o Wallac Oy Tel +358 2 267 8111 [maria.kivisto@pribori.pp.fi](mailto:maria.kivisto@pribori.pp.fi) **Scandinavia:** Reaction-Lab Oy Tel +358 2 410 1746 [info@reaction-lab.com](mailto:info@reaction-lab.com) **Spain** Izasa s.a. Tel +34 91 6630500 [cmaside@izasa.es](mailto:cmaside@izasa.es) **United Kingdom and U.S.A** Bioscan Inc Tel + 1 202 338 0974 [sales@bioscan.com](mailto:sales@bioscan.com)