


# UNICRYO FREEZE DRYER MC -60°C

- Vessel capacity 2 or 4 litre
- 750 to 1500 ml freezing process per 24 hours
- Transparent freezer drum
- Variable drying platforms with heated surface
- Defrosting heater
- Small footprint
- Fuzzy logic



Technische Änderungen ohne Vorankündigung auch innerhalb der Angebotsphase möglich  
Specifications are subject to change without notice for further improvements on design.

<http://www.uniequip.com>

UniEquip Products have CE-Test Certificate

# UNI-EQUIP

## UNICRYO MC2L-60°C FREEZE DRYER




Hedgehog with perspex drum

Freeze Drying Glassware  
50ml / 150ml / 300ml



Shelf Rack, 3 level platform  
with heating elements



Power Supply

Dimensions and specifications are subject to change without notice for further improvements on design.  
Dimensions may vary dependent on different accessories to be used.

# UNI-EQUIP

D-82152 Martinsried/München · Fraunhoferstrasse 11  
Telefon: +49-89-8 57 52 00 · Telefax: +49-89-8 56 13 04  
e-mail: uniequip@t-online.de · Internet: www.uniequip.com