

FILTERED LAMPS

> Power and versatility

Our elegant and highly efficient lamps are provided in 254, 312, 365 nm or combined. Power ranges from 4-watt up to 100-watt tube. Our long life filter has unlimited life for 312 and 365 nm and 3000 hours for 254 nm. Our unique filter minimizes white light interference allowing you to easily detect weak fluorescence.

Key features

- Single or dual wavelength
- Long life filter and high UV output
- Ondulex® reflector & aluminium housing for increased durability
- Easy to handle; stand or holder to add versatility

Specifications

Models	Tubes (Watt)	Wavelength (nm)	Intensity at 15cm (μW/cm ²)
VL-100.C	1 x 100-watt	254	3 400
VL-215.L	2 x 15-watt	365	2 300
VL-215.C	2 x 15-watt	254	1 780
VL-215.M	2 x 15-watt	312	3 000
VL-215.LC	1 x 15-watt	365	1 350
	1 x 15-watt	254	1 100
VL-215.LM	1 x 15-watt	365	1 350
	1 x 15-watt	312	1 800
VL-215.MC	1 x 15-watt	312	1 800
	1 x 15-watt	254	930
VL-115.L	1 x 15-watt	365	1 100
VL-115.C	1 x 15-watt	254	1 000
VL-115.M	1 x 15-watt	312	1 000
VL-8.L	1 x 8 -watt	365	800
VL-8.C	1 x 8 -watt	254	820
VL-8.M	1 x 8 -watt	312	790
VL-8.LC	1 x 8 -watt	365	720
	1 x 8 -watt	254	520
VL-8.LM	1 x 8 -watt	365	720
	1 x 8 -watt	312	660
VL-8.MC	1 x 8 -watt	312	660
	1 x 8 -watt	254	520
VL-6.L	1 x 6-watt	365	700
VL-6.C	1 x 6-watt	254	710
VL-6.M	1 x 6-watt	312	680
VL-6.LC	1 x 6-watt	365	610
	1 x 6-watt	254	400
VL-6.LM	1 x 6-watt	365	610
	1 x 6-watt	312	580
VL-6.MC	1 x 6-watt	312	580
	1 x 6-watt	254	400
VL-4.L	1 x 4-watt	365	400
VL-4.C	1 x 4-watt	254	340
VL-4.LC	1 x 4-watt	365	350
	1 x 4-watt	254	265

